

An ephys-compliant lick detector based on measured change in capacitance.

## Features

- High Frequency (100 or 125 [KHz]), low current (200[nA]) excitation signal makes this device invisible to [Neuropixel Probes](#) used in electrophysiology recordings.
- Fast. < 1[ms] response time.
- Contact-based. Device triggers when mouse tongue contacts either the dispensing tube *or* dangling reward liquid.
- TTL output triggers when a lick is detected.
- Harp-protocol compliant (serial num: 0x0578). Also dispatches timestamped Harp messages when lick state change has changed.
- Fully supported in Bonsai with a dedicated [Bonsai package](#)

## Extra Features

- 6-20VDC input (2.1 x 5.5mm barrel jack, positive center)
- reverse-polarity protected
- isolated USB to prevent ground loops with the PC.
- Two frequency options: 100 and 125[KHz] signal frequency to eliminate crosstalk between two closely-spaced lick tubes.
- Two amplitude options:
  - The 0.02[Vpp] option with proper grounding is intended for Ephys recordings as it introduces negligible noise artifacts.
  - The 2[Vpp] option can be used outside of an Ephys context without a ground connection

## Ordering

You can now order V1.0.0 from PCBWay:


- [main board](#). Select **PCB + Assembly**
- [cable buffer board](#). Select **PCB + Assembly**

These printed circuit boards are made on-demand.

You will also need:

- [Pololu 6-Pin SH Cable, 25cm long](#)
- [Stainless Steel Lick Tube, 18 Gauge](#) (any conductive lick tube will work!)
- [3D Printed Enclosure Walls and Laser Cut Lids](#)
  - right click the component to export each part.

## Device Pinout



## Switch Settings

This device has two settings (frequency and amplitude) that are independent of each other. These settings can be set by external switches to designate the state of the board upon power-up.

### ⚠ WARNING

Settings can be overwritten through software commands over the Harp interface after the board powers up.

## Low-Amplitude (Ephys) Setting

This setting uses the low-amplitude option of the signal.

It is **necessary** to use this setting when interacting with a mouse in an electrophysiology context (i.e: probes inserted in its brain).

Note that this setting, just like Neupixels Probes, is highly susceptible to noise. Good wiring hygiene (see below) is necessary to eliminate false positives from external noise sources.


## High Amplitude (Behavior) Setting

Outside of an electrophysiology context, this setting can be used. It makes the lick detector more immune to most common sources of EMI (valves, etc) because the signal is so much larger in comparison to the noise. This signal can also be run over a much longer cable. (We've used up to 2m!)


## Ephys Wiring Diagram

There are 3 configurations that will produce valid lick detection readings:


Option A works well if your rig *does not* have large sources of 60Hz noise.



Option B works well if your rig *does* have large sources of 60Hz noise.



Option C works well if your rig *does* have large sources of 60Hz noise, but you cannot earth ground it.



## Warnings

### Electrical Setup

It is critical that (1) both the device and mouse under test are grounded to a common ground and (2) the rig is earth-grounded (i.e: any exposed metal components are connected to earth ground). Otherwise, the device may introduce noise on the Neuropixel probes or produce spurious licks from outside electromagnetic interference.

### Multiple-Lick-Detector Crosstalk

Multiple lick detectors can be used in the same setup provided that (1) they are grounded correctly (see *Electrical Setup*), (2) the lick tubes are electrically insulated relative to each other, and (3) the spacing between exposed lick tubes exceeds 4mm. If the spacing between exposed metal lick tubes is less than 4mm, then the antenna noise from one lick detector may affect the other. This phenomenon is called *crosstalk*, and it can appear in 2 ways.

1. Lick signals that appear on one detector may show up on both detectors. The system is useable in this case, but the data may need to be cleaned up in post-processing where the shorter "simultaneously" detected lick is thrown out.
2. The amplitudes of the two AC signals will constructively interfere with each other, causing the excitation signals to saturate to their maximum values. The system cannot be used to detect licks this way.

To fix this issue,

1. Put each lick detector on a different frequency (100KHz and 125KHz).
2. Space the lick detection tubes farther apart.
3. Reduce the length of exposed metal on the lick tube.

# Theory of Operation

This device detects a threshold change in capacitance. A 100KHz, 10mVpp AC sine wave is played on the tip of a conductive lick spout, (such as McMaster-Carr pn: [89875K271](#)), and the amplitude of the returned waveform is measured every period.

A mouse sharing a ground with this device presents itself as a series resistive and capacitive load. When the mouse's tongue touches the metal tube (or water on the tip of the metal tube), the AC signal passes through the mouse to ground, lowering the amplitude of the returned sine wave below its trigger value, which triggers a lick to be detected. This detection signal is bandpassed in hardware (4th order Butterworth) and the resulting detection signal is lowpassed and compared against a threshold value in firmware to produce a lick/no-lick external trigger output. Overall propagation time through the entire signal chain is less than 1 millisecond.

The detection signal is AC such that repeated contact with the lick tube does not slowly charge the mouse.

## Tuning Parameters

**Warning:** changing the values of these tuning parameters may produce a detection signal that exceeds the <1[ms] detection time guarantee. We suggest ensuring that your setup is electrically grounded correctly first before changing these values.




There are three parameters that can be adjusted related to tuning the lick detection sensitivity. While the starting values should be sufficient for almost all setups, it is (unlikely but) possible that these knobs may need to be adjusted for a particular setup. The four parameters are:


- "fast moving average" window size
  - the last  $N$  raw signal amplitude values are averaged and sent to the consensus filter.
  - This filter smooths out high frequency amplitude changes by averaging the last few values together.
  - constraint:  $N$  must be a power of 2 (i.e: 2, 4, 8, 16, ...).
- "consensus" filter window size
  - the last  $N$  values must all be below the *lick-detection start* threshold to trigger a detected lick.
  - This filter increases the system's ability to reject intermittent (most likely inductive) noise sources that shrink the signal amplitude for a small period of time and may appear as a lick. Noise sources include fans, valves, and other lick detection tubes placed closer than 4mm apart.
  - Increasing this value will improve the system's ability to reject outside inductive noise, but will slow down the detection time.
  - constraint:  $N$  must be a power of 2.
- lick-detection "start" threshold
  - a signal value below this threshold will trigger a detected lick (i.e: the output signal will be set to 5V.)

- constraint: "start" < "stop" threshold
- lick-detection "stop" threshold
  - a signal value above threshold will untrigger a detected lick (i.e: the output signal will be set to 0V.)
  - constraint: "start" < "stop" threshold

LicketySplit	
whoAml	1400
firmwareVersion	0.0
hardwareTargets	0.5

## Registers

name	address	type	length	access	description	range	interfaceTyp
<a href="#">LickState</a> 	32	U8		Event	Emits an event when the state of any lick detector changes. Value will be High when lick detected and Low otherwise.		<a href="#">LickChannels</a> 
<a href="#">Channel0TriggerThreshold</a> 	33	U8		Write	Threshold value to detect the lick. Values below this threshold will be considered a detected lick.		

name	address	type	length	access	description	range	interfaceTyp
<a href="#">Channel0UntriggerThreshold</a> 	34	U8		Write	Threshold value to release the lick detection state. Values above this threshold will untrigger a detected lick.		

Analog Front-End for an implanted 22KOhm thermistor.

## Features:

- 1Hz highpass filter
- 1KHz lowpass filter
- $V_{out} = 1.65 + 100 \cdot V_{in}$
- Common mode thermistor input filtering
- Common mode supply voltage filtering

## Wiring Diagram

TODO

## Additional Parts

These parts aren't part of the PCBA, but they are required for a complete sensing solution:

- 22kOhm thermistor PN: GAG22K7MCD419 [digikey link](#)
- Flexible 2-conductor 36AWG Shielded TAU cable PN: TAU-3607.02 [link](#)
- 2-Pin Molex Connector Plug PN: 0533984002 ([link](#))
- 2-Pin Molex Connector Socket PN: 0151340206 ([link](#))

## Old Connectors

These connectors are no longer used:


- 2-pin Mill Max Connector pins PN: 850-10-002-20-001000 [digikey link](#)
- 2-socket Mill Max Connector socket PN: 801-47-002-10-012000 [digikey link](#)

SniffDetector	
whoAml	1401
firmwareVersion	0.1
hardwareTargets	0.1


## Registers

name	address	type	length	access	description	range	interfaceType
<a href="#">RawVoltage</a>	32	U16		Event	Emits a periodic event		



name	address	type	length	access	description	range	interfaceType
					containing the raw voltage read of the thermistor sensor.		
<a href="#">RawVoltageDispatchRate</a> 	33	U16		Write	Sets the rate at which the RawVoltage event is emitted.		

## A Harp Device for controlling an instrumented variable-friction mouse treadmill



### **NOTE**

This repository contains design files for the electronics assembly. For the mechanical assembly, design files, and bill-of-materials, see the [OnShape project](#).

## Features

- linear (no switching!) current control of 12V magnetic brake, Placid Industries B5ZD, [datasheet](#).
- reads measurements from an [RTS-10 torque transducer](#)
- wheatstone bridge input for torque transducer
- quadrature encoder input for position/speed measurements
- Harp-protocol compliant over full-speed USB
- accepts Harp time synchronization input.
- reverse-polarity protected power supply.

## Register Map

[List of Registers](#)





## Credits

This board design was created in partnership with the team at [Second Order Effects](#).

Treadmill	
whoAml	1402
firmwareVersion	0.0
hardwareTargets	0.0

## Registers

name	address	type	length	access	description	range	intei
<a href="#">Encoder</a>	32	S32		Read	Contains the current accumulated number of ticks.		
<a href="#">Torque</a>	33	S16		Read	Contains the current torque value.		
<a href="#">TorqueLoadCurrent</a>	34	S16		Read	Contains the current output current applied to the variable torque load.		
<a href="#">SensorData</a>	35	S32	3	Event	Emits a periodic event containing the packaged treadmill data. [Encoder, Torque, TorqueLoadCurrent]		<a href="#">Sens</a>
<a href="#">SensorDataDispatchRate</a>	36	U16		Write	Value greater than 0 will enable the periodic dispatch of treadmill data events at the specified rate (sp/s).	[0:1000]	
<a href="#">BrakeCurrentSetPoint</a>	37	U16		Write	Sets the raw value of the torque set-point to be applied to the	[0:65535]	

name	address	type	length	access	description	range	inter
					treadmill. This value is cleared to 0 if torque_limiting is enabled and triggered. Further writes in this condition return a WRITE_ERROR.		
<a href="#">TareSensors</a> 	38	U8		Write	Tares the specified sensors.		<a href="#">Sens</a>
<a href="#">ResetTareSensors</a> 	39	U8		Write	Removes the tare from the specified sensors.		<a href="#">Sens</a>
<a href="#">EnableTorqueLimit</a> 	40	U8		Write	Enables(1)/Disables(0) the brake if the maximum torque sensor value is detected. This register will be enabled by default.		<a href="#">Enab</a>
<a href="#">TorqueLimitState</a> 	41	U8		Event, Write	A value greater than 1 indicates that the torque limit has been triggered and the brake setpoint will be cleared. Writing a value of 0 will clear the torque limit state and re-enable the brake.		

An RP2040-based 16-output [Harp](#)-compliant Clock Synchronizer.

This device can serve as a drop-in replacement for the [Harp-Tech Clock Synchronizer](#) with added capabilities.



## Features

- 16 output channels for distributing clocks to other devices
- 1 input channel for receiving and synchronizing to another clock source
- 1 Additional Auxiliary Output, providing one of the following at a time:
  - PPS output
  - serial output of the time at user-specifiable baud-rate
- "Used channel detection." Device can identify which channels are in use.

## Ordering

You can now order V1.0.0 from PCBWay:

- [white-rabbit board](#). Select **PCB + Assembly**

These printed circuit boards are made on-demand.

You will also need:


- [3D Printed Enclosure Walls and Laser Cut Lids](#)
  - right click the component to export each part.

# Auxiliary Output

This device features an auxiliary output that can either produce a *pulse-per-second* (PPS) or UART message at the start of the whole second. This external signal enables Harp devices to additionally further synchronize with *non*-Harp devices.

## PPS Output

This device optionally outputs a 1[Hz] signal with a 50% duty cycle on the whole second.




Error from the nominal Harp time is < 1[us].


This feature is available on the AUX Port (3-pin terminal block).

## AUX UART Output

This device optionally outputs a the current time in seconds at the start of the whole second.



Error from the nominal Harp time is < 3[us].



This feature is available on the AUX Port (3-pin terminal block), and the baud rate is configurable via Harp Protocol (U32 in Register 36).

## PCBA Enclosure

For the enclosure design, see the companion [OnShape project](#).

## Developer Notes

In the firmware's **CMakeLists.txt**, adding (or uncommenting) the debug message below:

```
add_definitions(-DDEBUG)
```

will *override* the Auxiliary port behavior to use the auxiliary uart for `printf`-style debug messages. **PPS** and **AUX UART** features will not be available in debug mode.

WhiteRabbit	
whoAml	1404
firmwareVersion	0.1
hardwareTargets	0.1

## Registers

name	address	type	length	access	description	range	inter
<a href="#">ConnectedDevices</a>	32	U16		Event	The currently connected output channels. An event will be generated when any of the channels are connected or disconnected.		<a href="#">Clock</a>
<a href="#">Counter</a>	33	U32		Write, Event	The counter value. This value is incremented at the frequency specified by CounterFrequencyHz.		

name	address	type	length	access	description	range	inter
					Write to force a counter value.		
<a href="#">CounterFrequencyHz</a> <a href="#">↗</a>	34	U16		Write	The frequency at which the counter is incremented. A value of 0 disables the counter.	0 [0:500]	
<a href="#">AuxPortMode</a> <a href="#">↗</a>	35	U8		Write	The function of the auxiliary port.		<a href="#">AuxP</a> <a href="#">↗</a>
<a href="#">AuxPortBaudRate</a> <a href="#">↗</a>	36	U32		Write	The baud rate, in bps, of the auxiliary port when in HarpClock mode.	1000 [40:1000000]	



An octal TTL input-or-output board with programmable square wave output

## Features

- 8x TTL output or input, configurable.
- Programmable square wave generation.
- Harp-protocol compliant (serial num: 0x057B).
- Bonus: "passthrough buffer mode." External 3.3V and 5V CMOS devices can use this device as an octal buffer with external pins.

## Extra Features

- 6-20VDC input (2.1 x 5.5mm barrel jack, positive center)
- isolated USB to prevent forming ground loops with the PC.

## Wiring Diagram

TODO

## Theory of Operation

TODO: awesome diagram








## Harp Register Map









TODO

Cuttlefish	
whoAml	1403
firmwareVersion	0.1
hardwareTargets	0.0

## Registers

name	address	type	length	access	description	range	interfaceType
<a href="#">PortDirection</a> ↗	32	U8		Write	Set the direction of the ports		<a href="#">Ports</a> ↗

name	address	type	length	access	description	range	interfaceType
<a href="#">PortState</a> 	33	U8		Event, Write	Read or write the state of the ports. An event will be triggered when the state changes without a write command.		<a href="#">Ports</a> 
<a href="#">PwmTask</a> 	34	U8	18	Write	Struct to configure the PWM task. offset_us (U32), start_time_us (U32), stop_time_us (U32), port_mask (U8), cycles (U32),invert (U8)		
<a href="#">ArmExternalStartTrigger</a> 	35	U8		Write	If set to 1, the device will execute the PMW task using the selected pins.		<a href="#">Ports</a> 
<a href="#">ExternalStartTriggerEdge</a> 	36	U8		Write	Set the edge of the external		<a href="#">Ports</a> 

name	address	type	length	access	description	range	interfaceType
					trigger. 0: Rising, 1: Falling		
<a href="#">ArmExternalStopTrigger</a> 	37	U8		Write	If set to 1, the device will stop the PMW task using the selected pins.		<a href="#">Ports</a> 
<a href="#">ExternalStopTriggerEdge</a> 	38	U8		Write	Set the edge of the external trigger. 0: Rising, 1: Falling		<a href="#">Ports</a> 
<a href="#">SoftwareStartTrigger</a> 	39	U8		Write	Writing a non-0 value to this register will trigger the PWM task.		
<a href="#">SoftwareStopTrigger</a> 	40	U8		Write	Writing a non-0 value to this register will stop the PWM task.		
<a href="#">TaskControl</a> 	41	U8		Write			<a href="#">TaskControlPaylo</a> 

A self-contained, harp-compliant environment sensor

## Features

- pressure sensing
- temperature sensing
- humidity sensing
- USB powered
- Harp synchronization input

## Wiring Diagram

TODO

## Device Pinout

TODO

## Wiring Diagram

TODO

## Theory of Operation

TODO: awesome diagram







## Harp Register Map

[List of Registers](#)

EnvironmentSensor	
whoAml	1405
firmwareVersion	0.1
hardwareTargets	0.1

## Registers


name	address	type	length	access	description	range	interfaceType
<a href="#">Pressure</a>	32	U32		Read	Pressure, in Pa		
<a href="#">Temperature</a>	33	Float		Read	Temperature in		

name	address	type	length	access	description	range	interfaceType
					degrees C		
<a href="#">Humidity</a> 	34	Float		Read	Humidity, in %RH		
<a href="#">SensorData</a> 	35	Float	3	Event	A periodic event will be emitted with aggregated data from all sensors.		<a href="#">SensorDataPayload</a> 
<a href="#">EnableEvents</a> 	36	U8		Write	Enables (~2Hz) or disables the SensorData events		<a href="#">Events</a> 
<a href="#">TemperatureOffsetC</a> 	37	Float		Read	Fixed experimentally-determined calibration offset added to nominal temperature sensor reading in degrees C		

A [Cuttlefish](#) adapted to a dedicated fiber photometry application.

# Waveform Parameter Overview

This firmware generates the waveform below.



Multiple parameters are configurable:

- Continuous, or iterations
- software-trigger to start waveform generation
- **EVENT** message emitted upon rising edge of any camera trigger (configurable)


CuttlefishFip	
whoAml	1407
firmwareVersion	0.1
hardwareTargets	0.0

## Registers

name	address	type	length	access	description	range	interfaceType
<a href="#">SetTaskState</a>	32	U8		Write	Controls the state of tasks in the device.		<a href="#">TaskState</a>
<a href="#">AddTask</a>	33	U8	34	Write	Schedules a task by modelling following structure: U32 IOPin, float DutyCycle(0-1), float Frequency(Hz), U32 OutputMask (IOPins), U8 Events (0/1), U8 Mute (Kill the output but preserves timing), u32 delta1-4 (us)		
<a href="#">RemoveTask</a>	34	U8		Write	Removes a task scheduled at index 0-7. If the task is not scheduled, an error will be returned		<a href="#">TaskIndex</a>
<a href="#">ClearAllTasks</a>	35	U8		Write	Clears all scheduled task if a value of 1 is written.		<a href="#">EnableFlag</a>
<a href="#">TaskCount</a>	36	U8		Read	Returns the number of tasks currently scheduled. This		

name	address	type	length	access	description	range	interfaceType
					register is read-only.		
<a href="#">TaskRisingEdgeEvent</a> <a href="#">↗</a>	37	U8		Event	An event raised when a rising edge of any of the ports is detected. The <b>Events</b> flag must be enabled in the corresponding task to trigger this event. The event is raised when the task is started. The event is cleared when the task is removed or stopped.		<a href="#">Ports</a> <a href="#">↗</a>
<a href="#">Task0Settings</a> <a href="#">↗</a>	38	U8	34	Write	Represents the settings of Task0.		
<a href="#">Task1Settings</a> <a href="#">↗</a>	39	U8	34	Write	Represents the settings of Task1.		
<a href="#">Task2Settings</a> <a href="#">↗</a>	40	U8	34	Write	Represents the settings of Task2.		
<a href="#">Task3Settings</a> <a href="#">↗</a>	41	U8	34	Write	Represents the settings of Task3.		
<a href="#">Task4Settings</a> <a href="#">↗</a>	42	U8	34	Write	Represents the settings of Task4.		
<a href="#">Task5Settings</a> <a href="#">↗</a>	43	U8	34	Write	Represents the settings of Task5.		
<a href="#">Task6Settings</a> <a href="#">↗</a>	44	U8	34	Write	Represents the settings of Task6.		



name	address	type	length	access	description	range	interfaceType
<a href="#">Task7Settings</a> 	45	U8	34	Write	Represents the settings of Task7.		